**Caution:**

Don't mark on the phantom with pen or leave printed materials contacted on its surface.  
Ink marks on the phantom will be irremovable.

PH-50

# Newborn Whole Body Phantom

## Instruction Manual

### Contents

Please read	
General information .....	P.1
Training session	
Features • Training skills .....	P.2

*The world's first full body phantom for neonatal radiography.*

Newborn Whole Body Phantom is the world's first full body phantom for neonatal radiography with correct anatomical structure and movable limbs. Neonatal radiography is an important tool in NICU (Neonatal Intensive Care Unit). Patient positioning and immobilization are essential features. This phantom provides opportunities for hands-on training and experiments to minimize radiation exposure to newborn babies.

**Set includes**

Before your first use, please ensure that you have all components listed below.

Set includes: **a** newborn whole body phantom  
**b** storage case

1 **c** set of sample X-ray images 1  
1 instruction manual

Materials: Artificial bone: Epoxy resin  
Soft tissue: Urethane based resin

Phantom size:  
42 cm (represent a baby of 50 cm tall)/16.5 in  
Phantom weight: 2.8 kg /6.2 lb

Specifications are subject to change.

Anatomies -Skull, spine, clavicles, scapulae, ribs, pelvis, humeri, radius, ulnae, bones of hands, femora, fibulae, tibiae and bones of feet  
- Lungs, mediastinum, large intestine

**⚠ DOs and DON'T s**

**DO s**

Handle with care

It may cause the breakage in the phantom if the joints are forced to bend or degrees of angles outside of designed movable range.

Keep the phantom at room temperature, away from heat, moisture and direct sunlight.

**DON'Ts**

Never wipe the phantom or models with thinner or organic solvent.


Don't mark on the phantom with pen or leave printed materials contacted on their surface. Ink marks on the phantom will be irremovable.

Please note: The color of the phantom may change over time, though, please be assured that this is not deterioration of the material and the radiographic of the phantom stay unaffected.

### 1 Features • Training skills

#### ● Features

- Life size whole body newborn baby for plane radiography.
- Limbs rotate 360 degrees at shoulders and hip joints.
- Left hand is clenched and right hand is open.
- Original human tissue substitute.
- No metal parts or liquid structures.

#### ● Training skills

- Immobilization
  - Manual immobilization
  - Immobilization with fixtures
- Radiography
  - Upright PA/AP
  - Supine AP
  - Upright lateral
  - Supine lateral
- Autopsy imaging

Upright PA

Supine AP

## Caution

Don't mark on the phantom with pen or leave printed materials contacted on its surface.  
Ink marks on the phantom will be irremovable.

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