



Positioning of the body



Adjustable opening of the anus



Colonoscopy trainer for safe maneuvering

The soft, flexible colon tube provides a life-like response to the use of a colonoscope and facilitates training in skills such as preventing "loops" and safely advancing to the cecum.

Features

1. From simple and uncomplicated cases to difficult and sensitive variations. This training model is designed to allow for a wide range of difficulty for students and professionals alike.
2. The colon tube can be made airtight, allowing for air insufflation and suction.
3. The anal sphincter opening can be manipulated with a hand ball pump.

Training Skills

- Colonoscopy



Select difficulty level of training

Specifications :

Set Includes:		Manikin Size:	Replacement Parts:
1 abdomen model	1 plastic bottle	W48 x D22 x H30 cm	(11361-010) 1 colon-rectum tube
1 colon rectum tube	1 stand		(11361-020) 1 anus unit
1 anus unit	1 rubber sheet		(11361-030) 1 colon fixture set
1 air bulb	1 anus plug		(11361-120) 1 set of colon layout guide
1 abdomen skin cover	1 instruction manual		(11361-040) lubricant for colon-rectum tube (blue packs) 30 pcs x 50 ml each
1 set colon fixture	1 training guide		(11361-050) endoscope lubricant (pink packs) 30 pcs x 20 ml each
5 colon layout guides (6 cases total)	1 storage case		
30 lubricants for colon-rectum tube (50 ml each)			
30 endoscope lubricants (20 ml each)			

* Endoscope is not included.