



Two training pads to ensure improved skills in venipuncture

When working with Simulator Intravenous Arm II, trainees are to locate the vein by palpation to determine the puncture site and to recognize the needle tip's piercing the vein wall by sensation. The simulator allows for an array of quality skills training including venipuncture, IV access and injection.

Features

1. 2 kinds of puncture pads for difficulty variation
2. 3 accessible venous lines at cubital fossa:
Basilic vein, median vein, and cepahlic vein
3. Realistic needle tip resistance of the vein walls
4. "Flash-back" can be observed.

Training Skills

- Intravenous injection, blood collection and drip infusion



Specifications :

Set Includes:	Manikin Size:
2 arm models with supporter stand	W54 x D11 x H11 cm each
2 injection pads	
2 injection pads (red)	Replacement Parts:
1 circulation pump	(11267-300-03) 2 replacement injection pad for M50B
1 jar of coloring powder (red)	
1 set of spoon and plastic jar	(11267-300-04) 2 replacement injection pad for M50B (red)
1 instruction manual	