



### Product Supervision

Toyoaki Yamauchi, MD, ND, PhD  
 Professor  
 Department of Nursing Graduate School of Medicine Nagoya University




### Physiko can do what SP can not do.

Physiko is a physical assessment training simulator for medical and nursing students. Twelve preset examples of outpatients are designed to cover typical complaints, facilitating training in assessment procedures, assessment skills and providing hands-on experience of important abnormal findings. New patients for your original scenario can be created by simply selecting contents from built-in data.

### Cases

Assessment item	Condition		
Pupillary reflex	Normal	Pupillary constriction	
	Pupillary dilation	Bilateral asymmetry	
Blood pressure	Arbitrary-setup		
Auscultation of lung sounds	Normal	Bronchial breathing	Wheezes
	Weak in the left lung	Coarse crackles	Rhonchi
	Absent in the right lung	Fine crackles	
Auscultation of heart sounds	S2 split(-)	Aortic stenosis	Atrial fibrillation
	S2 split(+)	Mitral regurgitation	Arterial flutter
	S3 gallop	Mitral stenosis	Premature ventricular contraction (single)
	S4 gallop	Aortic regurgitation	Ventricular flutter
	S3 and S4 gallop	Sinus tachycardia	Ventricular fibrillation
	Innocent murmur	Sinus bradycardia	Sound volume adjuster
Auscultation of bowel sounds	Normal	Decrease	Ileus
	Increase	Subileus	
ECG simulation	Normal	Ventricular flutter	
	Atrial fibrillation	Ventricular fibrillation	
	Atrial flutter	Myocardial infarction (acute stage)	
	Premature ventricular contraction	Myocardial infarction (subacute stage)	
	Ventricular tachycardia	Myocardial infarction (chronic stage)	



\* A stethoscope and penlight are not included.


Patient Cases

Physiko's twelve preset patients cover typical outpatient complaints such as chest pain or breathing difficulties.

Top menu



Choose a patient



Start a session



Cases Comprehensive examination of 12 typical patients


Case	Age, Sex	Condition
Thoracic Pain	58, Male	Myocardial Infarction
	55, Male	Dissecting Aortic Aneurysm
	76, Female	Lung Infarction
	28, Male	Intercostal Muscle Ache
Abdominal Pain	26, Female	Ileus
	19, Male	Diarrhea

Case	Age, Sex	Condition
Lie Unconscious	36, Female	Brain Hypertension
Stuffy Feeling	84, Female	Pneumonia
	66, Male	Chronic Obstructive Lung Disease
	70, Male	Pulmonic Fibrosis
	73, Male	Heart Failure
	21, Female	Anemia

Physical Assessment Mode 2

Patient Cases Editor


Users can create a new simulated patient by selecting each findings from the built-in data.



Physical Examination

Skills Training

Physiko allows individual examination skills training.



Specifications :

Set Includes:			Manikin Size:
1 manikin	1 terminal box	1 ECG electrodes unit	155 cm
1 laptop PC	1 sphygmomanometer unit	1 set of clothes	