



### Score your suturing skills

Simple interrupted skin suturing can be evaluated quantitatively by validated evaluation items. Objective evaluation score and record of learning curve would encourage trainees to improve their skills effectively.

#### Features

1. There are three modes for different training purposes.
  - Evaluation: objective evaluation for self-learning
  - Learning: reference of examples and recorded data
  - Examination: skills assessment by the instructor
2. Six evaluation items leads trainees to desirable suture skills.

#### Training Skills

- Suture: Simple interrupted suture
- Ligature: Instrument tying

#### Evaluation Items

The sensors measure force and tension to the skin.

- Time
- Force on the tissue
- Suture tension


USB camera for Image Evaluation

- Spacing of stitches
- Equidistance
- Wound dehiscence

#### Specifications :

Set Includes:	Replacement Parts:
1 unit base	(11230-010)
1 skin suturing unit	40 simulated skins
40 sheets of simulated skin	
1 USB unit	
1 PC with software	
1 movie demonstration (CD)	
1 storage case	
1 instruction manual	

\* Instruments are not included.



### Quantitative assessment of suturing skills using laparoscopic forceps.

\*Training Box is not included in the standard set

#### Features

1. There are three modes for different training purposes.
  - Evaluation: objective evaluation for self-learning
  - Learning: reference of examples and recorded data
  - Examination: skills assessment by the instructor
2. Six evaluation items leads trainees to desirable suture skills.

#### Training Skills

- Suture: Simple interrupted suture
- Ligature: Instrument tying



#### Evaluation Items

The sensors measure force and tension to the skin.

- Time
- Force on the tissue
- Suture tension

USB camera for Image Evaluation

- Spacing of stitches
- Equidistance
- Wound dehiscence

#### Specifications :

Set Includes:
1 M57C suturing unit
40 sheets of simulated skin
1 USB unit
1 PC with software
1 demonstration movie (CD)
1 instruction manual
Replacement Parts:
(11230-010) 40 simulated skins

\* Instruments are not included.