

Strap-on Intramuscular Injection Trainer



The Strap-on Intramuscular Injection Trainer is designed for the practice between students and simulated patients. This trainer facilitates their understanding of the intramuscular injection process and the puncture feeling.

Features:

1. Durable sensor is embedded to evaluate the procedure.
2. Both lamp and buzz show success or failure in the procedure.
3. The water tank is located inside this trainer. It's compact and convenient.
4. This trainer is easier to drain without the waste water bag.
5. The skin and acromion is similar to the feeling of a real human and it allows palpating the anatomical landmark.



Strap-on Intramuscular Injection Trainer **M67B**

Set includes:

- | | |
|--------------------|--------------------|
| 1 Base Flam | 1 Waste Water Tank |
| 1 Skin | 1 Belt |
| 1 Soft Tissue Part | 1 Indication Box |
| 1 Muscular Part | 1 Power Adapter |
| 1 Acromion Part | |

Specifications are subject to change.

Replacement Parts:

- 11297-210
Replacement muscular part with sensor
- 11297-220
Replacement skin
- 11297-230
Replacement soft tissue part



11297-210
Replacement muscular parts with sensor



11297-220
Replacement skin



11297-230
Replacement soft tissue parts