

Breast Self Examination Model

Body-worn type



This model is effective and useful tool for Breast Self-Examination education. Newly developed material provides extremely human-like sensation.


The model facilitates practice in finding and distinguishing 6 different lumps through visual-tactile examination. It also offers good opportunities for training of professionals.

Features:

1. The model is designed to facilitate the diagnosis of six signs of breast tumors: skin dimples, cysts, cystoma, fibroadenoma, cancerous tumors.
2. This is perfect tool for breast self-examination education.
3. The appearance and texture of this life-size model are made as life-like as possible to facilitate clinical training.
4. The skin is easily replaceable.



A:mastopathy
B:cancer with depression in the skin
C:fibroadenoma
D:tumor
E:tumor
F:cancer with dimple sign



 **KYOTO KAGAKU co.,LTD**