



**Various types of skins to improve the basic of suturing skills**

**Features**

1. The artificial skin can be cut easily at any site to prepare a novel wound.
2. The soft skins are made of durable materials.
3. Simple and durable supporters keep the soft skins arched or maintain tension in the tubular models.
4. All instruments for the practice session are included in the package.

**Training Skills**

- 4 skills for suture: abdominal skin suture, blood vessel suture, intestinal suture, limbs skin suture



Abdomen



Blood vessel



Intestine



Limbs

**Specifications :**

Set Includes:		Replacement Parts:
1 set of simulated skins, ducts and supporters	1 set of instruments	(11219-030) 1 abdominal skin
- abdomen set: 2 abdominal skins, 1 supporter	- needle-holders	(11219-050) 5 blood vessels
- blood vessels: 5 narrow tubes	- suture needles	(11219-070) 5 intestinal tracts
- intestine set: 5 wide tubes, 1 supporter	- scissors	(11219-010) 5 limbs skins
- limbs set: 2 narrow skins, 1 supporter	- suture threads	
	- tweezers	
	1 storage case	