



Airway suction training for caregivers and healthcare professionals

A representation of human respiratory organs, designed for temporary airway suction and catheter insertion practice. Effective tool with life-like tactile qualities for training in airway suction.

Features

1. Realistic anatomy includes the mouth, tongue, oral and nasal pharynx, epiglottis, trachea and esophagus.
2. Catheterization pathways include oral, nasal or through the tracheostomy opening.
3. Life-size organs facilitate learning proper catheter insertion length.
4. The side window allows the passage of the catheter to be observed.
5. Realistic simulated sputum is sticky to the touch.
6. The trainer is washable and can be easily dismantled for cleaning.

Training Skills

- Temporary airway suction (intraoral and intranasal and endotracheal)

Nasal Suction



Oral Suction



Tracheal Suction



Specifications :

Set Includes:		Replacement Parts:	Manikin Size:
1 head & neck model with airway anatomy	1 metal tray	(11229-020) 1 replacement skin with trachea	W45 x D23 x H15 cm
	1 syringe	(11229-070) 1 simulated sputum	
1 bottle of lubricant	1 storage case	(11229-040) 2 sample catheters	
2 catheters	1 instruction manual	(11229-050) 1 lubricant	
1 bottle of simulated sputum			