



Quality and economy in one

Basic IV Arm is realistic with fine detail, inexpensive and durable, perfect for large number of trainees to practice skills required in their courses. Lightweight and portable, easy set up so that each trainee can have their own arm to practice. The vein tube can easily be filled with simulated blood using a soft plastic injection bottle.

Features

1. Flexible skin and palpable veins (replicates venous dilatation)
2. Realistic needle tip sensation and flashback
3. Injection pad is easily replaceable and durable.
4. No needle marks remain.
5. Hemorrhage can be stopped by compressing with the finger.
6. Adhesive plaster sticks well to the arm and can be easily removed.
7. Arm and supporter stand integrated so that users can set up easily.
8. Compact, lightweight, and portable

Training Skills

- Tourniquet application
- Venipuncture (including puncture with IV cannula)

Specifications :

Set Includes:	Manikin Size:
1 arm model with supporter stand	W54 x D11 x D11 cm
1 injection pad	
1 syringe (50ml)	
1 jar of coloring powder (red)	
1 spoon	
1 injection bottle for simulated blood	
1 carrying bag	



Affordable training unit for venipuncture

Strap-on Venipuncture Simulator is a simple, economical training kit for intravenous injection. This trainer is to be strapped onto the upper arm of an SP (Standardized Patient) or a life-size training manikin.

Features

1. Realistic tissue and vein wall resistance
2. Vein tube can be filled and refilled easily using a soft plastic bottle that works as a manual pump.
3. Three vein tubes are embedded: 2 normal and 1 thin.

Training Skills

- Blood collection and intravenous injection



Specifications :

Set Includes:
5 puncture pads (vein tube embedded)
5 supporting frames with belt
5 drain tubes
5 injector tubes
5 injection bottles
*Simulated blood powder is not included in the set.
Manikin Size:
W25 x D11 x H3 cm
Replacement Parts:
(11326-010) 1 replacement puncture pad
(11326-020) 5 replacement puncture pads

*Needle not included.