



Set of pediatric hands for practicing IV

Higher skills are essential in performing safe IV procedures on small infant hands with thin and soft veins that are difficult to find.

Features

1. The fingers and the back of the hand are flexible enough to be gripped or repositioned as in real-life procedures.
2. Two puncture sites on the back of the hand
3. Needle fixing procedure can be practiced.

Training Skills

- Pediatric dorsal veins injection
- Blood collection
- Intravenous infusion



Specifications :

Set Includes:
2 injection arm-hands (one-year-old)
2 injection arm-hands (three-year-old)
2 stands
1 circulation pump
1 jar of coloring powder (red)
1 spoon
1 plastic jar
1 plastic beaker
1 storage case
1 instruction manual
Manikin Size:
W23 x D6 x H5 cm (one-year-old)
W30 x D8 x H6 cm (three-year-old)
Replacement Parts:
(11350-010) 1 replacement injection arm-hand (three-year-old)
(11350-020) 1 replacement injection arm-hand (one-year-old)

* Needle not included.



Training in intradermal injection of forearm

This simulator provides training for intradermal injection to the inner forearm. Swelling of the skin will form when proper injection procedures are followed.

Features

1. Life-like touch and resistance of the injection site resembles the human arm
2. Injection pads are easily changeable.
3. No needle trace remains.

Training Skills

- Intradermal Injection



Specifications :

Set Includes:
5 puncture pads
1 forearm model
1 supporter stand
1 storage case
1 instruction manual
Manikin Size:
W20 x D20 x H8 cm
Replacement Parts:
(11346-010) 1 replacement puncture pad

* Needle not included.