



**Otoscope training with an alarm for painful insertions**

Ear Examination Simulator II is an upgraded trainer to acquire safe and reliable skills for examining the external acoustic meatus and tympanic membrane with a direct otoscope. Featured is an alarm when painful insertion is detected. The alarm has two modes: self-learning mode (with a beep), and instructor mode (silent).

**Features**

1. Alarm against painful insertion:
  - beep (self learning mode)
  - red light/yellow light (instructor mode)



**Training Skills**

- Ear examination with an otoscope
- Ear wax and foreign body removal

**Cases**

- Normal
- Serious otitis media (SOM)
- Mucoid otitis media (MOM)
- Chronic otitis media with perforation
- Acute suppurative otitis media (AOM)
- Cholesteatoma
- Tympanosclerosis
- Traumatic perforations
- Cerumen block



2. 9 major cases are preset.
3. Quick switch between cases.



4. Ear canals come in 2 sizes.
5. Anatomically correct soft ear model that allows training in straightening the auditory canal by pulling the auricle.
6. Work with a real otoscope.
7. The manikin neck can be tilted.
8. Training is possible with both left and right ears.

**Specifications :**

<b>Set Includes:</b>		<b>Replacement Parts:</b>
1 manikin head and shoulder	1 set of foreign body removal kit -simulated ear wax, small toys, small sponge fragments and a spoon	(11391-010) ear units and foreign body set -4 types of ear units - 1 pair each ear exam(normal), ear exam(stenosis), FB ear(normal), FB(stenosis) -1 set of foreign body removal kit simulated ear wax, small toys, small sponge fragments and a spoon
1 types of ear units -1 pair each -exam ear (normal) -exam ear(stenosis) -FB ear (normal) -FB ear (stenosis)		
1 instruction manual		
<b>Manikin Size:</b>		
W42 x D22 x H38 cm		