

2-in-1 simulator for prostate and rectal exams

The Prostate and Rectal Simulator includes 9 types of prostates and 4 types of rectums that are easily changeable. The lower torso manikin can be positioned in 3 different ways.

Features

1. Nine types of prostate simulate different scenarios
2. Four rectum units
3. Training with multiple prostates and rectal cases is concurrent and time efficient
4. Lateral, supine and prone positioning
5. 3 revolving units with 9 prostates can easily be interchanged.

Training Skills

- Digital-examination of the prostate and rectum
- Insertion and the use of the anal-scope and the proctoscope


Cases

Nine types of Prostates



- #1 normal
- #2 enlargement 1
- #3 enlargement 2
- #4 prostatitis
- #5 carcinoma 1
- #6 carcinoma 2
- #7 carcinoma 3
- #8 carcinoma 4
- #9 carcinoma 5

Four Rectum Units



- #1 normal
- #2 small carcinoma
- #3 large carcinoma
- #4 polyp combines with small carcinoma

Specifications :

Set Includes:	Manikin Size:
1 lower torso manikin	W34 x D33 x H24 cm
3 prostate units (include 9 types of prostates)	
4 rectum units	
1 Vaseline	
1 body	
1 storage case	