




With its various patient settings and a range of training features, "Yaye" brings uninterrupted training flow in community care and in-hospital procedures.

Featuring more than 40 training skills, "Yaye" serves an array of trainees from students to experienced nurses.

Training Skills



<p>Feeding</p> <ul style="list-style-type: none"> ● Tube Feeding (oral/nasal) <p>Elimination</p> <ul style="list-style-type: none"> ● Enema ● Urinary catheterization ● Stoma care ● Bowel irrigation ● Bedpan practice 	<p>Exercise</p> <ul style="list-style-type: none"> ● Body positioning ● Passive exercise ● Wheelchair transfer <p>Hygiene</p> <ul style="list-style-type: none"> ● Oral care ● Bed bathing ● Partial bathing ● Change of clothes ● Perineal care ● Hair care ● Grooming and appearance 	<p>Respiratory Care</p> <ul style="list-style-type: none"> ● Oxygen inhalation ● Airway suction ● Postural drainage ● Tracheotomy care <p>Wound Care</p> <ul style="list-style-type: none"> ● Pressure sore care ● Bandaging 	<p>Drug Administration</p> <ul style="list-style-type: none"> ● TPN care ● Suppository ● Subcutaneous injection ● Intramuscular injection ● Intravenous injection and infusion <p>Critical Care</p> <ul style="list-style-type: none"> ● Intubation assistance ● Chest compression ● Airway management 	<p>Assessment</p> <ul style="list-style-type: none"> ● Assessment of chest and abdomen <p>Pain Management</p> <ul style="list-style-type: none"> ● Fomentation <p>Others</p> <ul style="list-style-type: none"> ● Stomach pumping ● Vaginal douche ● Postmortem care
---	--	--	--	--

Features

ASSESSMENT OF CHEST AND ABDOMEN

Soft and elastic abdomen allows deep palpation



- Thoracoabdominal inspection, palpation and percussion
- Anatomical landmarks and soft tissue offer realistic training



APPEARANCE VARIETY

Patient variations with different age, sex and personality



Variations in face masks, wigs and external genitals offers wide possibilities of simulation settings and realistic training.



Easy Transformation

* The elderly face mask and the gray wig are optional parts.

Specifications:

Set Includes:		Replacement Parts:	
1 manikin	1 set of simulated suppository	(11404-020) face mask (young)	(11251-330) simulated feces set
1 denture	1 set of simulated feces	(11404-040) wig (young)	(11251-320) simulated suppository set
1 wig (black)	1 lubricant (for airway suction etc.)	(11404-060) female genitalia (catheterization, enema)	(11251-040) 2 thigh injection pads
1 face mask (young)	1 lubricant (for catheterization etc.)	(11251-310) female genitalia (suppository insertion, digital disimpaction)	(11251-030) 2 shoulder injection pads
1 injection pad (median vein)	1 irrigation bag	(11404-080)	(11251-050) 2 hip IM pads
2 shoulder injection pads	1 drainage pump (urinary bladder)	6 valves for urinary catheterization	(11404-050) 4 injection pad (median vein)
2 thigh injection pads	1 drainage pump (rectum)	(11251-350) rectal valve	(11229-050) lubricant (for airway suction etc.)
2 hip IM pads	1 drainage hose		(11381-800) lubricant (for catheterization)
1 female genitalia (catheterization)	1 set of pajamas		
1 female genitalia (suppository insertion/digital disimpaction)	1 instruction manual		
Manikin Size:		Optional Parts:	
152cm, 17kg		(11404-070) male genitalia	
		(11404-010) face mask (elderly)	
		(11404-030) wig (elderly)	

* Instruments are not included.

INFUSION CARE AND MANAGEMENT

Considerate IV route training



- **Left median vein:**
IV injection, infusion, securing cannula with tapes
- **Left hand dorsal vein:**
Securing cannula with tapes, caring during infusion



INTUBATION ASSISTANCE

Airway management and CPR



Training of intubation assistance can be performed in setting of critical care such as perioperative period and in ICU.

- Preparation of the equipment
- Oral intubation (laryngoscope and video laryngoscope)
- Confirmation of successful ventilation by:
 - Observation of thoracic movement (lung expansion)
 - Auscultation of the chest
- Securing the tube in place with tape
- Chest compression