



High quality models for catheterization training

High quality training in urethral catheterization, perineal care and enema. Soft and anatomically correct genitals and realistic feeling of catheter insertion. Organ unit can be used independently without the lumbar torso, allowing training with SPs or full body manikins in clinical settings. The modular system allows easy set-up, cleaning and maintenance, saving time and cost for training providers.

Features

1. Soft genitals allow realistic cleansing and other procedures such as holding the penis perpendicular to body.
2. Successful catheter insertion is confirmed by urine (water) outflow.
3. Manual bladder compression is possible.
4. Enema training in lateral position

Training Skills

- Intermittent catheterization
- Indwelling catheterization
- Perineal care
- Enema
- Localization of urethral meatus
- Cleansing
- Catheter insertion and placement
- Withdrawal of catheter
- Manual bladder compression



Enema Training



Passing catheter can be observed through transparent bladder reservoir.



Manual Bladder Compression



Training with SP

MW2C 11381-199-B

Genital unit and accessories for MW2A is sold as MW2C.



MW2D 11381-299-B

Genital unit and accessories for MW2B is sold as MW2D.



Specifications (MW2A or MW2B) :

Set Includes:		Manikin Size:	Replacement Parts:
1 lower torso manikin	1 rectum tube	W34 x D33 x H24 cm	(11381-300) 6 replacement valves for MW2A/B
1 genital organ unit	3 irrigation drainage tubes		(11381-610) 1 replacement female skin
3 valves (one is incorporated in the assembled unit)	1 lubricant		(11381-700) 2 replacement female meatus
	1 storage case		(11381-410) 1 replacement male skin
1 irrigation bottle	1 instruction manual		(11381-510) 2 replacement male meatus

* Catheters are not included.