

Recommended Devices

- Feeding tube: 5 fr.
- Suction catheter: 5 fr.
- Tracheal tube: 3.0mm ID
- Peripheral venous catheter: 28G x 30 cm
- Umbilical catheter: 5fr., 30cm
- Laryngoscope: #0
- BVM: 000



PI catheter insertion (vein localization using a transilluminator is possible)

NCPR Simulator Plus facilitates training of a variety of clinical skills required in NICU. The simulator is also suitable for the neonatal CPR training course.

Features

1. Training features of PI catheter insertion from right leg and arm are added to MW32 NCPR Simulator

Training Skills

- In addition to NCPR Simulator, NCPR Simulator Plus trains for:
 - ▶ IV access with transilluminator guidance
Basilic vein, dorsal vein (right arm/hand), Saphenous vein, popliteal vein (right leg)

Specifications :

Set Includes:
1 newborn manikin
1 lubricant
1 instruction manual

* Instruments are not included.



Recommended Devices

- Feeding tube: 3fr.
- Suction catheter: 3fr.
- Endotracheal tube: 2.0-2.5 mm ID
- Laryngoscope: #00
- Urine catheter: 6fr.
- BVM: neonatal (small)

Training Skills

- Changing body position, hygiene care in an incubator
- Assessment and measurement of whole body
- OG and NG Tube placement
- Insertion of urinary catheter
- BVM ventilation
- Chest compression
- Endotracheal intubation
- Nasal CPAP

*The product is not designed for bathing.

ELBW Simulator offers training opportunity for fundamental skills in caring an extremely low birth weight infant.

Features


1. A male infant of 26-29 weeks gestational, 900g weight
2. Resuscitation training including chest compression and endotracheal intubation
3. Incorporated anatomies such as bronchi, esophagus and stomach allows for training in basic skills
4. Articulated head and neck



Endotracheal intubation (confirmation of unilateral/bilateral ventilation)



Chest Compression



Procedures of tube feeding, oral and nasal airway suctioning

* Instruments are not included.

Specifications :

Set Includes:	Manikin Size:	Replacement Parts:
1 ELBW manikin	35 cm, 0.9kg	(11413-010) 1 mask for MW33
1 lubricant	Head circumference 28cm	(11229-050) 1 lubricant
1 instruction manual	Materials	
	Soft resin	