



Elastic external genitalia and realistic shape of the vaginal canal provides for a highly realistic training session, allowing the flow of skills from the perineal protection to the delivery of the placenta.

Features


1. Genital area made of innovative material that excels in elasticity, durability, and resilience, offers repetitive training.
2. Anatomically correct vaginal canal, ischial spine and pubic joint that facilitate training with anatomical understanding
3. Elastic eternal genitalia, realistic shape of the vaginal canal and the holding handle for the maternal body facilitate simulating a smooth delivery.

Training Skills

- Delivery assistance for different delivery positions (dorsal, all-fours)
- Perineal protection
- Delivery of the fetus (forceps delivery, vacuum extraction)
- Clamping, tying and cutting of the umbilical cords
- Delivery of the placenta
- Inspection of the velamen

Key Features

Realistic delivery of a baby with an elastic and resilient external genitalia



Skills training associated with anatomical understanding

Position of the fetus against the pelvis can be observed under direct vision during the delivery.



Two possible birth positions

- Dorsal position
- All-fours position

Specifications :

Set Includes:		Manikin Size:	Replacement Parts:	
1 maternal body	1 talcum powder	W74 x D44 x H33.5 cm 11.2kg	(11415-030) 1 genitalia unit	(11416-010) 1 fetus
1 genitalia unit (for pelvic examination and delivery assistance)	5 velamen sheets		(11416-020) 1 placenta	(11416-030) 10 umbilical cords for omphalotomy
1 fetus	3 lubricant (150 ml)	Materials	(11416-040) 5 velamen sheets	(11415-020) 3 lubricant (150 ml)
1 placenta	1 tray	Soft resin, hard plastic	Replacement Parts:	
3 umbilical cords for omphalotomy	1 instruction manual		(11417-020) 1 genitalia unit (perineorrhaphy) *skin included	