



Represents a perineal tear that occurs in a vaginal delivery. The innovative tissue material allows repetitive training.




Features

1. Suture training for first-degree perineal tears
2. Genital area that excels in elasticity, durability and resilience for multiple training opportunities.
3. Optional genitalia unit for training with the maternal body.

Training Skills

- Perineorrhaphy (perineal repair)
 - ▶ Single interrupted suture
 - ▶ Vertical mattress suture



First-degree perineal tear



* Suturing instruments are not included.



Genitalia skin can be mounted on the genitalia unit (optional part) for realistic training with the maternal body.

Specifications :

Set Includes:	Manikin Size:	Replacement Parts:
1 genitalia insert	W16 x D23 x H20.5 cm	(11417-010) 2 genitalia inserts (perineorrhaphy)
1 supporting holder	1.5kg	Optional Parts:
1 talcum powder	Materials	(11417-020) Genitalia unit (frame, skin)
1 instruction manual	Soft resin	