



Simulator for acquiring skills for epidural anesthesia that entail potential risks of complications

Epidural Anesthesia Simulator offers multiple training methods on both thoracic and lumbar areas, puncture block includes spinal columns from T7 to L5.

Features

1. Anatomically correct spinal column covering from T7 to L5.
2. Landmarks for palpation: IC, angulus inferior scapula and spinous process
3. True-to-life resistance and “pop” felt through the needles.
4. Left and right lateral positions
5. Loss of resistance technique is possible.
6. The broad puncture pad allows paramedian access as well as median approach.

Training Skills

- Epidural anesthesia procedures at lumbar and thoracic area
- Lumbar puncture




Landmark palpation



Epidural anesthesia



Lumbar puncture (CSF collection)



Specifications :

Set Includes:	Manikin Size:
1 life size unisex torso	W48 x D23 x H29 cm
5 puncture pads	Replacement Parts:
1 supporting stand	(11382-010)
1 50ml syringe	2 replacement puncture pads

* Needles are not included.