


## PH-10 BMD Chart Phantom UHA

*Bone Mineral Density chart for microdensitometry (MD) method.*



### Features


1. 21 steps with different hydroxyapatite content.
2. Steps range from 0 to 400 mg/cm<sup>2</sup>, with 20mg/cm<sup>2</sup> difference each.

### Specifications

Set Includes:	Size:
1 chart phantom	3 x 21 x 1.5 cm (consists of 21 blocks of 3 x 1 x 1.5 cm each.)
1 storage case	1.2 x 8.3 x 0.6 in (consists of 21 blocks of 8.3 x 0.4 x 0.6 in each)

## PH-16 Contrast Detail Phantom

*Image evaluation in plain X-ray*



### Features

- Four types of phantoms with different sizes and target types.
- Rod 15: 15 x 15 rods of height range from 1.0 to 8.0 mm (0.4 to 3.1 in)
  - Hole 15: 15 x 15 holes of depth range from 1.0 to 8.0 mm (0.4 to 3.1 in)
  - Rod 10: 10 x 10 rods of height range from 1.0 to 5.5 mm (0.4 to 2.2 in)
  - Hole 10: 10 x 10 holes of depth range from 1.0 to 5.5 mm (0.4 to 2.2 in)

### Specifications

Size:
Rod 15 and Hole 15: 24.5 x 24.5 cm each (9.6 x 9.6 in each)
Rod 10 and Hole 10: 17 x 17 cm each (6.7 x 6.7 in each)

## PH-17 Water Body Phantom WAC

*Water Body Phantom represents human chest and abdomen to serve as radiation absorber and scatterer.*




### Confirming to JIS Z 4915

### Specifications

Set Includes:	Size:
1 body phantom	30 x 20 x 45 cm
1 storage case	11.8 x 7.9 x 17.7 in

## PH-14 Acrylic Phantom XAC

*Slab phantoms for radiation absorption and scattering measurement.*



XAC-01	41430-000	30 x 30 x 0.1 cm/11.8 x 11.8 x 0.04 in
XAC-02	41431-000	30 x 30 x 0.2 cm/11.8 x 11.8 x 0.08 in
XAC-03	41432-000	30 x 30 x 0.3 cm/11.8 x 11.8 x 0.12 in
XAC-04	41433-000	30 x 30 x 0.4 cm/11.8 x 11.8 x 0.16 in
XAC-05	41434-000	30 x 30 x 0.5 cm/11.8 x 11.8 x 0.2 in
XAC-08	41435-000	30 x 30 x 0.8 cm/11.8 x 11.8 x 0.3 in
XAC-1	41436-000	30 x 30 x 1 cm/11.8 x 11.8 x 0.4 in
XAC-2	41437-000	30 x 30 x 2 cm/11.8 x 11.8 x 0.8 in
XAC-3	41438-000	30 x 30 x 3 cm/11.8 x 11.8 x 1.2 in
XAC-4	41439-000	30 x 30 x 4 cm/11.8 x 11.8 x 1.6 in
XAC-5	41440-000	30 x 30 x 5 cm/11.8 x 11.8 x 2.0 in
XAC-8	41441-000	30 x 30 x 8 cm/11.8 x 11.8 x 3.1 in
XAC-10	41442-000	30 x 30 x 10 cm/11.8 x 11.8 x 3.9 in