

PH-31

MRI Quality Assurance Phantom MHR

This QA phantom for MRI allows to evaluate the slice thickness, spatial resolution, uniformity and geometric distortion as well as contrast. Complies with NEMA standards.



Specifications

Set Includes:
1 phantom unit A
1 phantom unit B
1 set of nickel chloride solution

Size:
22 dia. x 14cm / 8.7 dia. x 5.5 in

PH-32

MRI Quality Assurance Phantom JMR

SN ratio, uniformity, slice thickness, spatial resolution, geometric distortion can be evaluated.



Conforming to JIS Z 4924

Specifications

Set Includes:
1 phantom
1 nickel chloride solution
1 storage case

Size:
22 dia. x 13 cm / 8.7 dia. x 5.1 in

PH-33

MRI Head Phantom NH

Life-size head phantom to assess uniformity.



Conforming to JIS Z 4924


Specifications

Set Includes:
1 head phantom
1 nickel chloride solution
1 storage case

PH-34

MRI/NM Head Phantom BHC

Simulate life-size head images in nuclear medicine and MRI.



Specifications

Set Includes:
1 head phantom
2 simulated tumor
(1 cm dia., 2 cm dia. each)
(0.4 in dia., 0.79 in dia. each)
1 nickel chloride solution
1 storage case

Size:
33 cm height
12.9 in height

Brain Phantom IB-20 advanced

This brain phantom of the striatal region with replicated skull densities of a male and female is useful for uptake ratio calibrations and studying the I-123 DaTSCAN scatter correction techniques.



Specifications

Size:
inside dimensions
13.5 x 18.5 cm / 5.3 x 7.3 in
height
8.8 cm / 3.5 in

