

## PH-40/41/42


# Tough Phantom Series

A stable, high quality and shatter-free phantom for radiotherapy.

## PH-40

### Tough Water Phantom WD

WD-3002	300 x 300 x 2 mm 12 x 12 x 0.08 in
WD-3003	300 x 300 x 3 mm 12 x 12 x 0.12 in
WD-3005	300 x 300 x 5 mm 12 x 12 x 0.2 in
WD-3010	300 x 300 x 10 mm 12 x 12 x 0.4 in
WD-3015	300 x 300 x 15 mm 12 x 12 x 0.6 in
WD-3020	300 x 300 x 20 mm 12 x 12 x 0.8 in
WD-3025	300 x 300 x 25 mm 12 x 12 x 1.0 in
WD-3030	300 x 300 x 30 mm 12 x 12 x 1.2 in
WD-3040	300 x 300 x 40 mm 12 x 12 x 1.6 in
WD-3050	300 x 300 x 50 mm 12 x 12 x 2.0 in
WD-4002	400 x 400 x 2 mm 16 x 16 x 0.08 in



WD-4003	400 x 400 x 3 mm 16 x 16 x 0.12 in
WD-4005	400 x 400 x 5 mm 16 x 16 x 0.2 in
WD-4010	400 x 400 x 10 mm 16 x 16 x 0.4 in
WD-4015	400 x 400 x 15 mm 16 x 16 x 0.6 in
WD-4020	400 x 400 x 20 mm 16 x 16 x 0.8 in
WD-4025	400 x 400 x 25 mm 16 x 16 x 1.0 in
WD-4030	400 x 400 x 30 mm 16 x 16 x 1.2 in
WD-4040	400 x 400 x 40 mm 16 x 16 x 1.6 in
WD-4050	400 x 400 x 50 mm 16 x 16 x 2.0 in

## PH-41

### Tough Bone Phantom BE-T, BE-H, BE-N

BE-T-2005	Compact Bone 200 x 200 x 5 mm / 8 x 8 x 0.2 in
BE-T-2010	Compact Bone 200 x 200 x 10 mm / 8 x 8 x 0.4 in
BE-T-2020	Compact Bone 200 x 200 x 20 mm / 8 x 8 x 0.8 in
BE-H-2005	Cortical Bone 200 x 200 x 5 mm / 8 x 8 x 0.2 in
BE-H-2010	Cortical Bone 200 x 200 x 10 mm / 8 x 8 x 0.4 in
BE-H-2020	Cortical Bone 200 x 200 x 20 mm / 8 x 8 x 0.8 in
BE-N-2005	Inner Bone 200 x 200 x 5 mm / 8 x 8 x 0.2 in
BE-N-2010	Inner Bone 200 x 200 x 10 mm / 8 x 8 x 0.4 in
BE-N-2020	Inner Bone 200 x 200 x 20 mm / 8 x 8 x 0.8 in



BE-H-3005	Cortical Bone 300 x 300 x 5 mm / 12 x 12 x 0.2 in
BE-H-3010	Cortical Bone 300 x 300 x 10 mm / 12 x 12 x 0.4 in
BE-H-3020	Cortical Bone 300 x 300 x 20 mm / 12 x 12 x 0.8 in
BE-N-3005	Inner Bone 300 x 300 x 5 mm / 12 x 12 x 0.2 in
BE-N-3010	Inner Bone 300 x 300 x 10 mm / 12 x 12 x 0.4 in
BE-N-3020	Inner Bone 300 x 300 x 20 mm / 12 x 12 x 0.8 in

## PH-42


### Tough Lung Phantom LP

LP-3010	300 x 300 x 10 mm / 12 x 12 x 0.4 in
LP-3020	300 x 300 x 20 mm / 12 x 12 x 0.8 in
LP-3030	300 x 300 x 30 mm / 12 x 12 x 1.2 in
LP-3050	300 x 300 x 50 mm / 12 x 12 x 2.0 in



### Easy-to-work

Tough series phantoms can be ordered with cavities and plugs.



## PH-37

# Therapy Body Phantom THRA-1

THRA-1 is an anthropomorphic, cross sectional dosimetry for therapeutic energy range.



### Specifications

Set Includes:  
1 torso phantom  
1 pair of breast phantom  
1 supporting frame  
1 storage case

Size:  
phantom size: 90 cm/35.5 in  
slice thickness: 3 cm/1.2 in

dosimeter holes:  
in lattice-like pattern of 3x3 cm/1.2 x 1.2 in

\*Slice thickness and dosimeter holes can also be custom ordered.

## PH-38

# Pediatric Therapy Body Phantom THRA-2



### Specifications

Set includes:  
1 torso phantom  
1 supporting frame  
1 storage case

Size:  
phantom size: 60 cm/23.6 in