



Thyroid Phantom

-Anthropomorphic Type-

PH-65



Anthropomorphic thyroid phantom simulates the absorption and scatter characteristics of human thyroid.

For quality assurance of system for iodine uptake ratio test; scatter, attenuation, and sensitivity.

Anatomies

- Cervical spine No.3~7
- Thoracic spine No.1
- Clavicle
- Breastbone (half)
- Thyroid (5 variations)

Variations of thyroid

	Volume
Thyroid 1	14.7 ml
Thyroid 2	16.7 ml
Thyroid 3	20.7 ml
Thyroid 4	30.2 ml
Thyroid 5	39.0 ml



Specifications

Set Includes:

- 1 neck and upper chest phantom
- 1 neck cover for thyroid
- 5 thyroid containers
- 1 storage case
- 1 instruction manual

Size:

- phantom size: 40 x 18 x 16 cm
- weight: 4 kg