

Ultrasound training for pediatric FAST procedures

This world's first pediatric ultrasound torso phantom provides opportunities of hands-on training in ultrasound examination of young patients.

Features


1. The phantom of a 2-year-old, thoracoabdominal organs, a bone structure, free-fluid to learn FAST procedures and pathologies that are commonly seen in pediatrics.
2. With this phantom, trainees can acquire skills in the basics of pediatric abdominal ultrasound as well as FAST/Acute sonography.

Training Skills

- FAST procedures:
Internal hemorrhage at perihepatic, perisplenic, pelvis, and pericardium area
- FAST procedures:
Internal hemorrhage at perihepatic, perisplenic, pelvis, and pericardium area

Product Supervision

Shunsuke Nosaka, MD, PhD
 Director of Radiology
 National Medical Center for Children and Mothers



Specifications :

Set Includes:	Phantom Size:
1 ultrasound phantom	W41 x D15 x H15 cm
1 storage case	
1 tutorial manual	