Prostate Examination Simulator


Prostate Examination Simulator provides realistic training in prostate palpation skills, with its soft life-like anus and rectum material.

Ten types of prostates are placed in the simulator and the trainer can smoothly interchange them using a simple rotating system.

Three patient positions are possible. Compact and space-saving design.

Production supervision:

Dr. Osamu Yoshida, Professor Kyoto University Faculty of Medicine, Department of Urology Dr. Takashi Katayama, Professor Toyama Medical and Pharmaceutical University Department of Urology













The lumbar torso can be placed in three patient positions: supine, prone and lateral position.

Built-in prostates can be changed using a simple rotating switch.

This saves time for preparation in training and skills assessment.



Prostate Examination Simulator

Set includes:

M 53B

- 1 male lumbar torso
- 1 rotating unit
 - 10 prostate cases: 2 normal, 2 enlargements, 5 carcinomas and 1 prostatitis
- 1 jar of Vaseline
- 1 pair of gloves
- 1 carrying case

manikin size: $40 \times 32 \times 24 \text{ H cm}$, 5.3 kg packing size: $41 \times 40 \times 33 \text{ H cm}$, 7 kg

Replacement parts

11272-040 replacement anus pad

Specifications are subject to change.

KYOTO KAGAKU

15 Kitanekoya-cho, Fushimi-ku, Kyoto, 612-8388, JAPAN Tel: +81-75-605-2510 Fax: +81-75-605-2519 www.kyotokagaku.com rw-kyoto@kyotokagaku.co.jp