

# Strap-on Venipuncture Trainer

Product No. **M87**

## Set includes:

- 5 puncture pads (vein tube embedded)
- 5 supporting frames with belt
- 5 injector tubes
- 5 drain tubes
- 5 injection bottles
- \*Simulated blood powder is not included in the set.

## Replacement parts

- 11326-010 replacement puncture pad (1 piece)
- 11326-020 replacement puncture pad (a set of 5)

Specifications are subject to change.

Simple, economical training kit for daily training.


Vein tube makes three lines, two normal and one thin, around the center of the pad.

Realistic tissue and vein wall resistance.

Vein tube can be filled and refilled easily using soft plastic bottle that works as a manual pump.

Puncture pad is replaceable.

# Intradermal Injection Training Simulator

Product No. **M94**

## Set includes:

- 5 puncture pads
  - 1 forearm model
  - 1 stand
  - 1 carrying case
- packing size: 31 x 34 x 14 H cm 3.3kg

## Replacement parts

- 11346-010 replacement puncture pad

Specifications are subject to change.

This simulator provides training of intradermal injection to the inner forearm.

A clear characteristic skin wheal is formed when proper injection procedure is followed.

Life-like touch and resistance of the injection site, close to that of a real human arm, allows pulling the skin and making punctures like a real patient.

Life-size forearm model made of soft resin, realistic in appearance and touch, help in training to find a proper needle angle.

Stand supports the forearm model at proper position.

Special structure of the injection pad, when properly injected, creates a clear characteristic skin wheal.

Injection pads are easily exchangeable.

No needle trace remains.

 **KYOTO KAGAKU**

www.kyotokagaku.com rw-kyoto@kyotokagaku.co.jp

Head Office: 15 Kitanechoya-cho, Fushimi-ku, Kyoto, 612-8388, JAPAN Tel: +81-75-605-2510 Fax: +81-75-605-2519

USA Office: 3109 Lomita Boulevard, Torrance, CA 90505-5108, USA Tel: 1-310-325-8860 Fax: 1-310-325-8867