

Digital Rectal Examination Simulator

Digital Rectal Examination Simulator is a unique trainer designed for training in differentiation of rectal tumors and polyps by palpation.


Four interchangeable rectum units are included.

Normal prostate or endocervix is also to be palpated.

Life-like rectum is palpable up to seven cm depth.

High-fidelity sphincter resistance and feeling of rectum wall allows excellent hands-on DRE training before passing on to patients.


4 interchangeable rectum units for differentiation of lesion types.
Prostate or endocervix can be palpated.



Prostate or endocervix can be palpated.



4 interchangeable rectum units for differentiation of lesion types.



Digital Rectal Examination Simulator

Set includes:

M 92

- 1 lumbar torso (unisex)
- 4 rectum units:
 - 1 normal
 - 2 rectal cancers
 - 1 rectum with polyps
- 1 prostate unit
- 1 endocervix unit
- 1 jar of Vaseline
- 1 carrying case

manikin size: 30 x 35 x 29 H cm, 6.2 kg
packing size: 44 x 53 x 53 H cm, 13 kg

Specifications are subject to change.

 KYOTO KAGAKU

15 Kitanechoya-cho, Fushimi-ku, Kyoto, 612-8388, JAPAN
Tel: +81-75-605-2510 Fax: +81-75-605-2519
www.kyotokagaku.com rw-kyoto@kyotokagaku.co.jp

Replacement parts

- 11269-110 rectum unit #1 (normal)
- 11269-120 rectum unit #2 (cancer A)
- 11269-130 rectum unit #3 (cancer B)
- 11269-140 rectum unit #4 (rectum with polyps)
- 11269-150 prostate unit
- 11269-160 endocervix unit