


Peripheral Venous Catheter Placement Simulator

MW9



median antebrachial vein



dorsal vein

Peripheral Venous Catheter Placement Simulator provides full procedural skills from injection to fixing the IV cannula.

Features

1. All training skills for the peripheral venous route management
2. Two puncture sites: median antebrachial vein and dorsal vein
3. Practice for the extension of hand skin and the angiopressure management
4. Durable puncture pad with IV cannula

Training Skills

- Wrapping with tourniquet
- Confirmation of puncture site
- Sterilization of puncture site
- Puncture with IV cannula
- Confirmation of back flow in puncture
- Angiopressure management and decannulation
- Setting of infusion tube
- Confirmation of natural instillation
- Fixing of the puncture site
- Injection of medical solution from injection sub port



Puncture with IV cannula




Confirmation of flash back in puncture



Angiopressure management and decannulation



Confirmation of natural instillation



Fixing of the puncture site



Injection of medical solution from injection sub port

Set Includes

- 1 arm model
- 2 injection pads for dorsal vein of hand
- 2 injection pads for median antebrachial vein
- 1 bottle for simulated blood
- 1 stand for bottle of simulated blood
- 1 infusion set
- 1 infusion bag
- 1 syringe (50ml)
- 1 plastic beaker
- 10 simulated blood sticks (swab type)
- 1 instruction manual

Manikin Size

W52 x D10 x H10 cm

Replacement Parts

- 11388-200
- 2 injection pads for median antebrachial vein
- 11388-300
- 2 injection pads for dorsal vein of hand
- 11388-400
- 10 simulated blood sticks (swab type)