

# Pitting Edema Model

Five variations of pitting edema models including a normal one.  
Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"  
The models can be worn on the lower thighs of a SP or a training manikin.

Production Supervision:  
Junji Otaki ProfessorTokyo Medical University



Four grades of pitting edemas with different resistances, depth of indentations.

## Degrees of pitting edemas:



**Normal**

No indentation left



**+ Light pitting**

Indentations with approx. 2mm depth.



**2+ Moderate pitting**

Indentations with approx. 4mm depth.



**3+ Deep pitting**

Veins, bones and joints are indistinct because of swell.  
Indentations with approx. 6mm depth.




**4+ Very deep pitting**

Indentations with approx. 8mm depth.

## Training items:

Differentiation of pitting edemas with four different grades and a normal.  
Learn the examination procedures for edemas on the lower thighs.  
The removable fibula plate facilitates training without SP leg or simulation on the soft bases.



## Pitting Edema Model

**MW 1**

Set includes:  
strap-on pitting edema model  
5 variations 1 piece each  
packing size: 26 x 23 x 12H cm, 1 kg  
Specifications are subject to change.

**KYOTO KAGAKU**

15 KitaneKoya-cho, Fushimi-ku Kyoto,  
612-8388, JAPAN  
Tel: +81-75-605-2510  
Fax: +81-75-605-2519  
www.kyotokagaku.com  
rw-kyoto@kyotokagaku.co.jp