

Pitting Edema Model

Five variations of pitting edema models including a normal one.
Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"
The models can be worn on the lower thighs of a SP or a training manikin.

Production Supervision:
Junji Otaki ProfessorTokyo Medical University



Four grades of pitting edemas with different resistances, depth of indentations.

Degrees of pitting edemas:



Normal

No indentation left



+ Light pitting

Indentations with approx. 2mm depth.



2+ Moderate pitting

Indentations with approx. 4mm depth.



3+ Deep pitting

Veins, bones and joints are indistinct because of swell.

Indentations with approx. 6mm depth.




4+ Very deep pitting

Indentations with approx. 8mm depth.

Training items:

Differentiation of pitting edemas with four different grades and a normal.
Learn the examination procedures for edemas on the lower thighs.

The removable fibula plate facilitates training without SP leg or simulation on the soft bases.



Pitting Edema Model

MW 1

Set includes:
strap-on pitting edema model
5 variations 1 piece each
packing size: 26 x 23 x 12H cm, 1 kg
Specifications are subject to change.

KYOTO KAGAKU

15 KitaneKoya-cho, Fushimi-ku Kyoto,
612-8388, JAPAN

Tel: +81-75-605-2510

Fax: +81-75-605-2519

www.kyotokagaku.com

rw-kyoto@kyotokagaku.co.jp