

Ultrasound Examination Training Phantom "ECHOZY"

Ultrasound Training Phantom with no pathology

Upper abdomen organs with approximate-to-human echogenicities.
 State-of-art anatomy fulfills medical sonographers' high requirements.
 Unique long-life phantom materials.



Production supervision:
 Dr. Hitoshi Asai, Director
 Osaka Kyoiku University, Health
 Administration Center
 Dr. Shigeru Nakamura
 Nagayoshi General Hospital, Clinical
 Examination Department



Training items

Instrumentation and probe manipulation

Basics of abdominal sonography:


- Cross sections and sonographic anatomy
- Sonographic demonstration of each individual organ
- Localization of hepatic Couinaud's segments

M 98

Internal Organ Anatomical Model "ECHO-ZOU"

Separates into 20 parts:

1. Hepatic segment
2. Hepatic segment
3. Hepatic segment
4. Hepatic segment
5. Hepatic segment
6. Hepatic segment
7. Hepatic segment
8. Hepatic segment
9. Portal vein, Bile duct and Hepatic vein
10. Gallbladder
11. Pancreas
12. Spleen
13. Right kidney
14. Left kidney
15. Abdominal aorta
16. Inferior vena cava
17. Hepatic veins
18. Spinal column
19. Stomach
20. Large intestine, Small intestine







□ KYOTO KAGAKU

Ultrasound Examination Training Model

Set variations:

*please quote the eight digit product number when you place orders.

US-1B	41900-030	"ABDFAN"
<p>1 ultrasound phantom "ABDFAN" 1 set of positioning pillows 1 tutorial Manual (DVD) 1 carrying case</p> <p>phantom size: approx.18 x 25 x 28H cm,12 kg packing size: approx.30 x 70 x 50H cm,17 kg</p>		
US-1B	41900-100	"ABDFAN" Full Set
<p>1 ultrasound phantom "ABDFAN" 1 anatomical model "ECHO-ZOU" 1 set of positioning pillows 1 tutorial Manual (DVD) 1 carrying case</p> <p>phantom size: approx.18 x 25 x 28H cm,12 kg ECHO-ZOU size: approx.16 x 23 x 16H cm, 2.9 kg packing size: approx.30 x 70 x 50H cm,17 kg</p>		
US-1	41900-010	"ECHOZY"
<p>1 ultrasound phantom "ECHOZY" 1 set of positioning pillows 1 carrying case</p> <p>phantom size: approx.18 x 25 x 28H cm,12 kg packing size: approx.30 x 70 x 50H cm,17 kg</p>		
US-1	41900-000	"ECHOZY" Full Set
<p>1 ultrasound phantom "ECHOZY" 1 anatomical model "ECHO-ZOU" 1 set of positioning pillows 1 carrying case</p> <p>phantom size: approx.18 x 25 x 28H cm,12 kg ECHO-ZOU size: approx.16 x 23 x 16H cm, 2.9 kg packing size: approx.30 x 70 x 50H cm,17 kg</p>		

Materials:
 ultrasound phantom: Urethane elastomer

 **KYOTO KAGAKU**

15 Kitanechoya-cho, Fushimi-ku, Kyoto, 612-8388, JAPAN
 Tel: +81-75-605-2510 Fax: +81-75-605-2519
 www.kyotokagaku.com rw-kyoto@kyotokagaku.co.jp