

M92



Digital Rectal Examination Simulator

This Digital Rectal Examination Simulator is a unique trainer designed for palpation training to identify rectal tumors and polyps.



FEATURES

- 1 | 4 interchangeable rectum units for identification of lesion types.
- 2 | Normal prostate or endocervix can also be palpated.
- 3 | Life like rectum is palpable up to 7cm deep
- 4 | High-fidelity sphincter resistance and feeling of rectum wall allow for excellent hands-on DRE training before examining actual patients.

SKILLS

- | Digital palpation of the prostate

DESCRIPTIONS

SPECIFICATIONS

Size: W37 x D28 x H30 cm
W14.5 x D11 x H11.8 inch
Weight: 6kg / 13.3lbs

SET INCLUDES

1 lumbar torso (unisex)	1 prostate unit
4 rectum units:	1 endocervix unit
1 normal	1 Vaseline
2 rectal cancers	1 storage case
1 rectum with polyps	1 instruction manual

REPLACEMENT PARTS

11269-210	1 rectum unit #1 (normal)
11269-220	1 rectum unit #2 (cancer A)
11269-230	1 rectum unit #3 (cancer B)
11269-240	1 rectum unit #4 (rectum with polyps)
11269-150	1 prostate unit
11269-160	1 endocervix unit

M53B



Prostate Examination Simulator

Prostate Examination Simulator provides realistic training in prostate palpation skills with its soft life-like anus and rectum material.



FEATURES

- 1 | Ten cases of prostate are set on a rotating table inside the body
- 2 | The cases are easily changeable
- 3 | Three patient positions are possible: supine, prone and lateral position

SKILLS

- | Digital palpation of the prostate

CASES

- | 2 normal, 2 enlargements, 5 carcinomas, 1 prostatitis

DESCRIPTIONS

SPECIFICATIONS

Size: W30 x D36 x H24 cm
W11.8 x D14.2 x H9.4 inch

Weight: 5.1kg/ 10.9lbs

MATERIALS

Soft resin, Hard resin
Latex free

SET INCLUDES

1 male lumbar torso
1 rotating unit
10 prostate cases: 2 normal,
2 enlargements, 5 carcinomas
and 1 prostatitis
1 lubricant gel
1 pair of gloves
1 storage case
1 instruction manual

REPLACEMENT PARTS

11272-040	1 anus pad
-----------	------------