


M85



Airway Suction Trainer

Airway suction and tracheostomy management simulator

Life-like anatomical and tactile features to further understanding temporary airway suctioning



Nasal suction



Oral suction

FEATURES

- 1 | Facilitate understanding of the way to insert catheter.
- 2 | Transparent cross-sectional window to monitor the catheter
- 3 | Realistically viscous simulated sputum is included
- 4 | Can be easily dismantled for cleaning
- 5 | Able to open the mouth by hand

SKILLS

- | Oral suctioning
- | Nasal suctioning
- | Endotracheal suctioning

LANDMARKS

- | oral cavity
- | tongue
- | oral and nasal pharynx
- | epiglottis
- | trachea
- | esophagus

OPTIONAL PARTS

11229-060 Oral Cover for training using liquid

Oral cover to store simulated sputum can be purchased as optional parts. This oral cover is to maintain the sputum inside the oral cavity.



DESCRIPTIONS

SPECIFICATIONS

Size: 45 x 23 x 15 cm / 17.7 x 9 x 5.9 inch

Weight: 2kg/4.4lbs

SET INCLUDES

1 head & neck model with airway anatomy	1 syringe
1 bottle of lubricant	1 storage case
2 catheters	1 lubricant
1 bottle of simulated sputum	1 instruction manual
1 metal tray	

MATERIALS

Soft resin

Latex free

REPLACEMENT PARTS

11229-020	1 replacement skin with trachea
11229-070	1 simulated sputum
11229-040	2 sample catheters
11229-050	1 lubricant

OPTIONAL PARTS

11229-060	oral cover for suctioning
-----------	---------------------------