

MW30

# Head and Neck Palpation Trainer

## New, innovative and ultra-realistic touch, designed to standardize palpation experiences

With realistic tactile material, this head and neck palpation trainer provides a sensation similar to human skin

### Product Supervision

Junji Otaki,  
professor and acting chief of the center,  
Center for Medical Education,  
Hokkaido University



### KEY FEATURES



#### Swallowing

Swallowing motion to identify and assess the thyroid gland



#### Sites

Various deviations



#### Close-to-real-life sensation

Supple, resilient and delicate texture

### FEATURES


- 1 | Provides opportunities for hands-on training.
- 2 | Abnormal findings to complement training with healthy simulated patients
- 3 | Ensures that trainees share the same tactile sensation to specify the findings
- 4 | Maintain objective approaches for skills examination

### SKILLS

- Palpation of head and neck for lymph nodes, thyroid and salivary glands
- Palpation skills to support clinical reasoning

### CASES

Sites	Findings
1   Posterior auricular nodes	Edema/Hard
2   Occipital nodes	N/A
3   Superficial cervical nodes	Edema/Palpable
4   Posterior cervical nodes	Edema/Palpable
5   Deep cervical nodes	Edema/Palpable
6   Pre-auricular nodes	Edema/Hard
7   Tonsillar nodes	Edema/Palpable
8   Submandibular nodes	Edema/Palpable
9   Submental nodes	Edema/Palpable
10   Supraclavicular nodes	Edema/Hard/Immovable
11   Subclavicular (infraclavicular) nodes	Edema/Hard/Immovable
12   Thyroid gland	Graves' disease



## MW31 Head and Neck Palpation Trainer for OSCE

The palpable nodes and glands in the standard type are distributed between two identical manikins for use in situations like student exams or skill checking

### Internal findings

- A** left side: 6/1/3/10  
right side: 4/5/13/12
- B** left side: 3/13/ 9/5  
right side: 4/11/5

### DESCRIPTIONS

#### SPECIFICATIONS

Size: W29 x D29 x H47 cm / W11.4 x D11.4 x H18.5 in. Weight: 6.5kg / 14.3lbs

#### SET INCLUDES

- 1 | head and neck manikin with storage case
- 1 | instruction manual

#### MATERIALS

Soft resin  
Latex free