





Anorectal Manometry Examination Trainer

Realistic model for training in anorectal transducer insertion



11460-000



11460-100

FEATURES

(common to 11460-000 & 11460-100)

- 1 | Designed specifically for insertion training of rectal pressure transducers
- 2 | The simulated pressure transducer with scale marking reproduces the flexibility of an actual device, providing a realistic insertion feel.

SKILLS

- | Insertion of rectal transducer
- | Confirmation of insertion depth

11460-000 Rectal Manometry Examination Model

- 1 | Life-sized buttocks that allows learning consideration about patient discomfort in addition to skills of transducer insertion
- 2 | The anal unit is detachable, allowing confirmation of the insertion depth

11460-100 Training Set of 5

- 1 | Compact and space-saving design, ideal for group training
- 2 | Transparent rectum allows visual confirmation of insertion depth
- 3 | Recommended to be used together with 11460-000 for comprehensive training including patient discomfort

DESCRIPTIONS

SET INCLUDES

11460-000

- 1 lumbar torso
- 1 rectal unit (pre-attached to main body)
- 1 simulated pressure transducer
- 1 petroleum jelly (Vaseline)
- 1 talcum powder
- 1 instruction manual

11460-100

- 5 table-top training unit
- 5 simulated pressure transducer
- 1 petroleum jelly (Vaseline)
- 1 talcum powder
- 1 instruction manual

SPECIFICATIONS

11460-000

Size: W31 × D28 × H30 cm
Weight: 2.8 kg

MATERIALS

Soft resin, hard plastic
Latex free

11460-100

Size: W13 × D17 × H12.5 cm (per model)
Weight: 0.7 kg × 5 (total 3.5 kg)

REPLACEMENT PARTS

11460-010 Simulated Pressure Transducer