



OASIS Perineorrhaphy Simulator -4th degree-

Develop diagnostic and restorative skills for OASIS (Obstetric Anal Sphincter Injuries)

Product Supervision

Dr. MAKI Jota, Assistant Professor, Department of Obstetrics and Gynecology, Graduate School of Biomedical Sciences, Okayama University




fig. Structure of the MW81


FEATURES

- 1 | The OASIS condition, in which the rectal mucosa has been lacerated down to the rectum, is reproduced. Mucosal and muscular layers are arranged in the same way as in the human body.
- 2 | Palpation and suture sensations are reproduced in the same way as in the human body.
- 3 | The genitalia skin is designed to be removable and can be easily replaced after wear and tear.

SKILLS

Perineorrhaphy

- | Rectovaginal Fistula
- | Rectal Mucosa
- | Internal Anal Sphincter
- | External Anal Sphincter
- | Vagina Mucosa



The genital skin is designed for easy attachment and removal.

DESCRIPTIONS

SPECIFICATIONS

Size: W160×D140×H180 mm

Weight: 0.8 kg

SET INCLUDES

OASIS simulator:	1	instruction manual
1 supporting holder		
1 genitalia skin		
1 external anal sphincter		
1 fixing clamp		

MATERIALS

Soft resin, Latex free

REPLACEMENT PARTS

11485-010 1 genitalia skin (OASIS 4th degree perineorrhaphy)