

MW82 Femoral Arterial Puncture Trainer

Hands-on femoral arterial access training both with ultrasound-guided and landmark-guided procedures.

Features

- Realistic arterial pulse
- Palpable ASIS and pubic bone landmarks
- Physiological counterflow
- No visible puncture marks at the injection site
- Leak-proof connections



Skills

- Ultrasound-guided femoral artery puncture
- Landmark-guided femoral artery puncture
- Blood sampling

Femoral Arterial Puncture Trainer

MW82

11486-000

Set Includes

groin manikin	1
circulation pump (AC100V-240V)	1
artery tube	2
vein tube	2
cotton swab red	10
instruction manual	

Specifications

	Groin manikin
dimensions	W33 x D35 x H15 cm (manikin) W20 x D16 x H13 cm (pump)
weight	3.4 kg (manikin) 1.5kg (pump)
materials	Soft resin, Latex free



Recommended devices

Needle size: 23G

Replacement parts



Replacement skin
11486-030



artery tube (a set of 4)
11486-010



vein tube (a set of 4)
11486-020



cotton swab red (set of 10)
11388-400

Related products



An essential simulator to increase confidence and competence

Arterial Puncture Training Wrist
M99 11351-000



Master Ultrasound-Guided Radial Arterial Puncture

Ultrasound Guided Arterial Puncture Training Pads
11351-110



A sequence of uninterrupted simulation experience for specialist training

CVC Insertion Simulator II & III
M93UB, U, C 11347-300, 500, 510

KYOTO KAGAKU

<https://www.kyotokagaku.com/>
rw-kyoto@kyotokagaku.co.jp

2026.03

Kyoto Head Office / Factory / Show Room

15 Kitanekeya-cho Fushimi-ku Kyoto, Japan 612-8388
Telephone:+81-75-605-2510 Fax:+81-75-605-2519

Asia and other regions

Overseas Div. / MEM Training Center

Second fl. NREG Hongo-3-chome bldg., 3-26-6 Hongo,
Bunkyo-ku, Tokyo 113-0033, Japan
Telephone:+81-3-3817-8071 Fax:+81-3-3817-8075

All American regions

Kyoto Kagaku America Inc. / Show Room

3109 Lomita Boulevard, Torrance, CA 90505-5108, USA
Telephone:+1-310-325-8860 Fax:+1-310-325-8867
Email : info@kkamerica-inc.com

Europe, Russia & Africa

Kyoto Kagaku Europe GmbH / Training Room

De-Saint-Exupery-Str.10, 60549 Frankfurt, Germany
Telephone:+49-69-5060-28160