

US-21

# DVT (Deep Vein Thrombosis) Phantom

For ultrasound screening of DVT (Deep Vein Thrombosis)

## Product Supervision

Kenji Karino, Professor, Clinical Skill-Up Center, Simane University Hospital





## FEATURES

- 1 | Faithfull reproduction of blood vessel meandering, branching, and thrombus.
- 2 | The built-in blood vessels allows for milking method as well as observation of compression of the vessels by applying pressure with a transducer.
- 3 | Familiarize with assessment of venous thrombosis in the lower limbs.

## SKILLS

- | Screening for deep venous thrombosis in the lower limbs
- | Observation of blood flow by milking method

## Images



## DESCRIPTIONS

### SPECIFICATIONS

Size: W27x D17 x 47H cm each  
W10.6x D6.7x H18.5inch each

Weight: 3.0kg each  
6.6lbs each

### SET INCLUDES

- 1 lower limb manikin
- 1 ultrasound calf phantom (incorporated)
- 1 talcum powder
- 1 storage case
- 1 instruction manual

\* Ultrasound scanner is not included.